



UBA1 Antibody, Biotin conjugated

Product Code	CSB-PA00149D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22314
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 1 protein (730-881AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like modifier-activating enzyme 1 (EC 6.2.1.45) (Protein A1S9) (Ubiquitin-activating enzyme E1), UBA1, A1S9T UBE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.