




AMPD1 Antibody

Product Code	CSB-PA001680LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23109
Immunogen	Recombinant Human AMP deaminase 1 protein (97-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AMP deaminase 1 (EC 3.5.4.6) (AMP deaminase isoform M) (Myoadenylate deaminase), AMPD1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	AMPD1

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, PC-3 whole cell lysate, A549 whole cell lysate

All lanes: AMPD1 antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 91, 90 kDa

Observed band size: 91 kDa