




APOA5 Antibody

Product Code	CSB-PA001917LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Q788
Immunogen	Recombinant Human Apolipoprotein A-V protein (260-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein A-V (Apo-AV) (ApoA-V) (Apolipoprotein A5) (Regeneration-associated protein 3), APOA5, RAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOA5

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA001917LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.