



ATP1B1 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA002326LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05026 |
| Immunogen | Recombinant Human Sodium/potassium-transporting ATPase subunit beta-1 protein (64-241AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sodium/potassium-transporting ATPase subunit beta-1 (Sodium/potassium-dependent ATPase subunit beta-1), ATP1B1, ATP1B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ATP1B1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |