



ATP1B1 Antibody, FITC conjugated

Product Code	CSB-PA002326LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05026
Immunogen	Recombinant Human Sodium/potassium-transporting ATPase subunit beta-1 protein (64-241AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/potassium-transporting ATPase subunit beta-1 (Sodium/potassium-dependent ATPase subunit beta-1), ATP1B1, ATP1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATP1B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.