



BMP10 Antibody, Biotin conjugated

Product Code	CSB-PA002734LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95393
Immunogen	Recombinant Human Bone morphogenetic protein 10 protein (317-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 10 (BMP-10), BMP10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BMP10