



C1S Antibody, Biotin conjugated

Product Code	CSB-PA003657LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09871
Immunogen	Recombinant Human Complement C1s subcomponent protein (519-688AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1s subcomponent (EC 3.4.21.42) (C1 esterase) (Complement component 1 subcomponent s) [Cleaved into: Complement C1s subcomponent heavy chain; Complement C1s subcomponent light chain], C1S
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	C1S
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.