



CCRL2 Antibody, FITC conjugated

Product Code	CSB-PA004852LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00421
Immunogen	Recombinant Human C-C chemokine receptor-like 2 protein (1-43AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor-like 2 (Chemokine receptor CCR11) (Chemokine receptor X) (Putative MCP-1 chemokine receptor), CCRL2, CCR11 CCR6 CKRX CRAM HCR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCRL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.