



CD3G Antibody, Biotin conjugated

Product Code	CSB-PA004933LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09693
Immunogen	Recombinant Human T-cell surface glycoprotein CD3 gamma chain protein (138-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD3 gamma chain (T-cell receptor T3 gamma chain) (CD antigen CD3g), CD3G, T3G
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD3G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.