



CKAP4 Antibody, Biotin conjugated

Product Code	CSB-PA005453LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07065
Immunogen	Recombinant Human Cytoskeleton-associated protein 4 protein (260-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoskeleton-associated protein 4 (63-kDa cytoskeleton-linking membrane protein) (Climp-63) (p63), CKAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	CKAP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.