



CMKLR1 Antibody, FITC conjugated

Product Code	CSB-PA005614LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99788
Immunogen	Recombinant Human Chemokine-like receptor 1 protein (321-373AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chemokine-like receptor 1 (G-protein coupled receptor ChemR23) (G-protein coupled receptor DEZ), CMKLR1, CHEMR23 DEZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CMKLR1