



CYP2A6 Antibody, HRP conjugated

Product Code	CSB-PA006410LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11509
Immunogen	Recombinant Human Cytochrome P450 2A6 protein (86-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 2A6 (EC 1.14.13.-) (1,4-cineole 2-exo-monooxygenase) (CYP2A6) (Coumarin 7-hydroxylase) (Cytochrome P450 2A6) (Cytochrome P450(I)), CYP2A6, CYP2A3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP2A6