



CYP2A6 Antibody, HRP conjugated

Product Code	CSB-PA006410LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11509
Immunogen	Recombinant Human Cytochrome P450 2A6 protein (86-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 2A6 (EC 1.14.13.-) (1,4-cineole 2-exo-monooxygenase) (CYP2A6) (Coumarin 7-hydroxylase) (Cytochrome P450 IIA3) (Cytochrome P450(I)), CYP2A6, CYP2A3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP2A6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.