





CYP2B6 Antibody

Product Code	CSB-PA006414LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20813
Immunogen	Recombinant Human Cytochrome P450 2B6 protein (250-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 2B6 (EC 1.14.13.-) (1,4-cineole 2-exo-monooxygenase) (CYP11B6) (Cytochrome P450 IIB1), CYP2B6
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP2B6

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA006414LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA006414LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.