




DHODH Antibody

Product Code	CSB-PA006852LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02127
Immunogen	Recombinant Human Dihydroorotate dehydrogenase (quinone), mitochondrial protein (241-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydroorotate dehydrogenase (quinone), mitochondrial (DHODEHase) (EC 1.3.5.2) (Dihydroorotate oxidase), DHODH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DHODH

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA006852LA01HU at dilution of 1:100