



DHODH Antibody, HRP conjugated

Product Code	CSB-PA006852LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02127
Immunogen	Recombinant Human Dihydroorotate dehydrogenase (quinone), mitochondrial protein (241-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydroorotate dehydrogenase (quinone), mitochondrial (DHODEHase) (EC 1.3.5.2) (Dihydroorotate oxidase), DHODH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DHODH