




DHPS Antibody

Product Code	CSB-PA006854LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49366
Immunogen	Recombinant Human Deoxyhypusine synthase protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyhypusine synthase (DHS) (EC 2.5.1.46), DHPS, DS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DHPS

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA006854LA01HU at dilution of 1:100