



DUSP1 Antibody, Biotin conjugated

Product Code	CSB-PA007238LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28562
Immunogen	Recombinant Human Dual specificity protein phosphatase 1 protein (125-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 1 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH1) (Mitogen-activated protein kinase phosphatase 1) (MAP kinase phosphatase 1) (MKP-1) (Protein-tyrosine phosphatase CL100), DUSP1, CL100 MKP1 PTPN10 VH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.