



E2F4 Antibody, Biotin conjugated

Product Code	CSB-PA00759D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16254
Immunogen	Recombinant Human Transcription factor E2F4 protein (160-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor E2F4 (E2F-4), E2F4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	E2F4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.