



EMX2 Antibody, Biotin conjugated

Product Code	CSB-PA007657LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04743
Immunogen	Recombinant Human Homeobox protein EMX2 protein (75-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein EMX2 (Empty spiracles homolog 2) (Empty spiracles-like protein 2), EMX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EMX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.