



F13A1 Antibody, HRP conjugated

Product Code	CSB-PA007919LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00488
Immunogen	Recombinant Human Coagulation factor XIII A chain protein (380-520AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor XIII A chain (Coagulation factor XIIIa) (EC 2.3.2.13) (Protein-glutamine gamma-glutamyltransferase A chain) (Transglutaminase A chain), F13A1, F13A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F13A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.