





FGF17 Antibody

Product Code	CSB-PA008622LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60258
Immunogen	Recombinant Human Fibroblast growth factor 17 protein (23-216AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 17 (FGF-17), FGF17
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGF17

Image



Western Blot
 Positive WB detected in: Mouse liver tissue
 All lanes: FGF17 antibody at 3µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 25, 24 kDa
 Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA008622LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.