



FGF22 Antibody, Biotin conjugated

Product Code	CSB-PA008628LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCT0
Immunogen	Recombinant Human Fibroblast growth factor 22 protein (23-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 22 (FGF-22), FGF22
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGF22
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.