



FLNA Antibody, Biotin conjugated

Product Code	CSB-PA008724LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21333
Immunogen	Recombinant Human Filamin-A protein (1234-1462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Filamin-A (FLN-A) (Actin-binding protein 280) (ABP-280) (Alpha-filamin) (Endothelial actin-binding protein) (Filamin-1) (Non-muscle filamin), FLNA, FLN FLN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FLNA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.