




FOXH1 Antibody

Product Code	CSB-PA008817LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75593
Immunogen	Recombinant Human Forkhead box protein H1 protein (153-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein H1 (Forkhead activin signal transducer 1) (Fast-1) (hFAST-1) (Forkhead activin signal transducer 2) (Fast-2), FOXH1, FAST1 FAST2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXH1

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, 293 whole cell lysate

All lanes: FOXH1 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 40 kDa

Observed band size: 40 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.