



GNPAT Antibody, HRP conjugated

Product Code	CSB-PA009629LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15228
Immunogen	Recombinant Human Dihydroxyacetone phosphate acyltransferase protein (589-680AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydroxyacetone phosphate acyltransferase (DAP-AT) (DHAP-AT) (EC 2.3.1.42) (Acyl-CoA:dihydroxyacetonephosphateacyltransferase) (Glycerone-phosphate O-acyltransferase), GNPAT, DAPAT DHAPAT
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GNPAT