



GRIK1 Antibody, HRP conjugated

Product Code	CSB-PA009906LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39086
Immunogen	Recombinant Human Glutamate receptor ionotropic, kainate 1 protein (675-834AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, kainate 1 (GluK1) (Excitatory amino acid receptor 3) (EAA3) (Glutamate receptor 5) (GluR-5) (GluR5), GRIK1, GLUR5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.