



HLTF Antibody, Biotin conjugated

Product Code	CSB-PA010520LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14527
Immunogen	Recombinant Human Helicase-like transcription factor protein (332-476AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Helicase-like transcription factor (EC 2.3.2.27) (EC 3.6.4.-) (DNA-binding protein/plasminogen activator inhibitor 1 regulator) (HIP116) (RING finger protein 80) (RING-type E3 ubiquitin transferase HLTF) (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 3) (Sucrose nonfermenting protein 2-like 3), HLTF, HIP116A RNF80 SMARCA3 SNF2L3 ZBU1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HLTF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.