





HSF2 Antibody

Product Code	CSB-PA010792LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03933
Immunogen	Recombinant Human Heat shock factor protein 2 protein (79-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock factor protein 2 (HSF 2) (Heat shock transcription factor 2) (HSTF 2), HSF2, HSTF2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HSF2


Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA010792LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA010792LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA010792LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)