



HTR7 Antibody, HRP conjugated

Product Code	CSB-PA010899LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34969
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 7 protein (1-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 7 (5-HT-7) (5-HT7) (5-HT-X) (Serotonin receptor 7), HTR7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HTR7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.