



# IGFBP6 Antibody, HRP conjugated

Product Code	CSB-PA011100LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24592
Immunogen	Recombinant Human Insulin-like growth factor-binding protein 6 protein (32-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-like growth factor-binding protein 6 (IBP-6) (IGF-binding protein 6) (IGFBP-6), IGFBP6, IBP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IGFBP6