



IRF9 Antibody, Biotin conjugated

Product Code	CSB-PA011824LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00978
Immunogen	Recombinant Human Interferon regulatory factor 9 protein (44-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 9 (IRF-9) (IFN-alpha-responsive transcription factor subunit) (ISGF3 p48 subunit) (Interferon-stimulated gene factor 3 gamma) (ISGF-3 gamma) (Transcriptional regulator ISGF3 subunit gamma), IRF9, ISGF3G
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.