



ITGA1 Antibody, FITC conjugated

Product Code	CSB-PA011861LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56199
Immunogen	Recombinant Human Integrin alpha-1 protein (242-524AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-1 (CD49 antigen-like family member A) (Laminin and collagen receptor) (VLA-1) (CD antigen CD49a), ITGA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.