



ITGA6 Antibody, HRP conjugated

Product Code	CSB-PA011869LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23229
Immunogen	Recombinant Human Integrin alpha-6 protein (625-866AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-6 (CD49 antigen-like family member F) (VLA-6) (CD antigen CD49f) [Cleaved into: Integrin alpha-6 heavy chain; Integrin alpha-6 light chain; Processed integrin alpha-6 (Alpha6p)], ITGA6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.