



ITGAV Antibody, Biotin conjugated

Product Code	CSB-PA011877LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06756
Immunogen	Recombinant Human Integrin alpha-V protein (775-885AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-V (Vitronectin receptor) (Vitronectin receptor subunit alpha) (CD antigen CD51) [Cleaved into: Integrin alpha-V heavy chain; Integrin alpha-V light chain], ITGAV, MSK8 VNRA VTNR
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ITGAV
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.