



KIF5B Antibody, Biotin conjugated

Product Code	CSB-PA012342LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33176
Immunogen	Recombinant Human Kinesin-1 heavy chain protein (338-514AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-1 heavy chain (Conventional kinesin heavy chain) (Ubiquitous kinesin heavy chain) (UKHC), KIF5B, KNS KNS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.