



KIF5B Antibody, Biotin conjugated

Product Code	CSB-PA012342LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33176
Immunogen	Recombinant Human Kinesin-1 heavy chain protein (338-514AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-1 heavy chain (Conventional kinesin heavy chain) (Ubiquitous kinesin heavy chain) (UKHC), KIF5B, KNS KNS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF5B