



LAMB1 Antibody, FITC conjugated

Product Code	CSB-PA012730LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07942
Immunogen	Recombinant Human Laminin subunit beta-1 protein (523-661AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit beta-1 (Laminin B1 chain) (Laminin-1 subunit beta) (Laminin-10 subunit beta) (Laminin-12 subunit beta) (Laminin-2 subunit beta) (Laminin-6 subunit beta) (Laminin-8 subunit beta), LAMB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LAMB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.