




LTA4H Antibody

Product Code	CSB-PA013219LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09960
Immunogen	Recombinant Human Leukotriene A-4 hydrolase protein (107-311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukotriene A-4 hydrolase (LTA-4 hydrolase) (EC 3.3.2.6) (Leukotriene A(4) hydrolase), LTA4H, LTA4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LTA4H

Image



Western Blot

Positive WB detected in: 293 whole cell lysate


All lanes: LTA4H antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 70, 60, 58, 67 kDa

Observed band size: 70 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA013219LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.