



LTA4H Antibody, FITC conjugated

Product Code	CSB-PA013219LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09960
Immunogen	Recombinant Human Leukotriene A-4 hydrolase protein (107-311AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukotriene A-4 hydrolase (LTA-4 hydrolase) (EC 3.3.2.6) (Leukotriene A(4) hydrolase), LTA4H, LTA4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LTA4H
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.