



MAF Antibody, HRP conjugated

Product Code	CSB-PA013312LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75444
Immunogen	Recombinant Human Transcription factor Maf protein (55-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor Maf (Proto-oncogene c-Maf) (V-maf musculoaponeurotic fibrosarcoma oncogene homolog), MAF
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.