



MAST3 Antibody, Biotin conjugated

Product Code	CSB-PA013513LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60307
Immunogen	Recombinant Human Microtubule-associated serine/threonine-protein kinase 3 protein (1101-1309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated serine/threonine-protein kinase 3 (EC 2.7.11.1), MAST3, KIAA0561
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAST3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.