



MCM3AP Antibody, Biotin conjugated

Product Code	CSB-PA013592LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60318
Immunogen	Recombinant Human Germinal-center associated nuclear protein (292-435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Germinal-center associated nuclear protein (GANP) (80 kDa MCM3-associated protein) (MCM3 acetylating protein) (MCM3AP) (EC 2.3.1.-) (MCM3 acetyltransferase), MCM3AP, GANP KIAA0572 MAP80
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MCM3AP