





MPP1 Antibody

Product Code	CSB-PA014758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00013
Immunogen	Recombinant Human 55 kDa erythrocyte membrane protein (123-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	55 kDa erythrocyte membrane protein (p55) (Membrane protein, palmitoylated 1), MPP1, DXS552E EMP55
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPP1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA014758LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA014758LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.