



MPP1 Antibody, FITC conjugated

Product Code	CSB-PA014758LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00013
Immunogen	Recombinant Human 55 kDa erythrocyte membrane protein (123-252AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	55 kDa erythrocyte membrane protein (p55) (Membrane protein, palmitoylated 1), MPP1, DXS552E EMP55
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPP1