



MPST Antibody, FITC conjugated

Product Code	CSB-PA014770LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25325
Immunogen	Recombinant Human 3-mercaptopyruvate sulfurtransferase protein (102-208AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-mercaptopyruvate sulfurtransferase (MST) (EC 2.8.1.2), MPST, TST2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPST
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.