





NCAM2 Antibody

Product Code	CSB-PA015512LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15394
Immunogen	Recombinant Human Neural cell adhesion molecule 2 protein (209-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neural cell adhesion molecule 2 (N-CAM-2) (NCAM-2), NCAM2, NCAM21
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NCAM2


Image



Western Blot
Positive WB detected in: MCF-7 whole cell lysate
All lanes: NCAM2 antibody at 2.7µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 94, 48 kDa
Observed band size: 94 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA015512LA01HU at dilution of 1:100



Immunofluorescence staining of HepG2 cells with CSB-PA015512LA01HU at 1:266, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.