



NCAM2 Antibody, FITC conjugated

Product Code	CSB-PA015512LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15394
Immunogen	Recombinant Human Neural cell adhesion molecule 2 protein (209-333AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neural cell adhesion molecule 2 (N-CAM-2) (NCAM-2), NCAM2, NCAM21
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NCAM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.