



NCR3 Antibody, Biotin conjugated

Product Code	CSB-PA015551LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14931
Immunogen	Recombinant Human Natural cytotoxicity triggering receptor 3 protein (19-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural cytotoxicity triggering receptor 3 (Activating natural killer receptor p30) (Natural killer cell p30-related protein) (NK-p30) (NKp30) (CD antigen CD337), NCR3, 1C7 LY117
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NCR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.