



NOX1 Antibody, Biotin conjugated

Product Code	CSB-PA015959LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5S8
Immunogen	Recombinant Human NADPH oxidase 1 protein (418-564AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH oxidase 1 (NOX-1) (EC 1.-.-.) (Mitogenic oxidase 1) (MOX-1) (NADH/NADPH mitogenic oxidase subunit P65-MOX) (NOH-1), NOX1, MOX1 NOH1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NOX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.